SYNOLIS VA HYALURONIC ACID 40mg / SORBITOL 80mg

GB Instructions for use of SYNOLIS VA 40/80 p
FR Notice d'utilisation de SYNOLIS VA 40/80 p

p2 p7

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Tel: +41 (0) 848 949 848 mail@aptissen.com www.aptissen.com

Geneva/Switzerland

Updated documentation is available from Aptissen SA. in Switzerland.

 (ϵ)

1639

PSL010-4

Description

Synolis VA is a viscoelastic, sterile, apyrogenic, isotonic, buffered, 2% solution of sodium hyaluronate. Sodium hyaluronate used in Synolis VA is obtained from bacterial fermentation and has a high mean molecular weight (m.w.) of 2 MDaltons. Performance of Synolis VA in the treatment of symptomatic osteoarthritis (OA) is achieved through high concentration and high mean molecular weight of sodium hyaluronate combined with an agent which limits HA degradation (sorbitol), when tested in vitro under non-physiological conditions. This unique combination confers to this viscoelastic solution its ability to restore joint lubrication and its shock absorbing properties similar to those of the healthy synovial fluid.

Synolis VA functions by restoring the physiological and viscoelastic properties of the synovial fluid which has progressively been lost during the OA development. herefore, Synolis VA reduces the local pain and discomfort caused by OA, and it improves mobility of the synovial joints.

Introduction

Synolis VA is presented in a pre-filled 2ml glass syringe. The syringe is packaged in a sealed blister. Each box contains 1 instruction leaflet.

Composition

For 1 mL:

Sodium hyaluronate 20mg Sorbitol 40mg Ruffered saline solution OSP 1ml

One syringe contains 2mL of viscoelastic solution. Sterilisation method: Moist heat

Indications

Synolis VA is indicated for the treatment of pain and functional impairment as a consequence of degenerative changes (osteoarthritis) of the knee.

Dosage and method of administration

Synolis VA should be injected within the synovial cavity by a physician skilled in performing intra-articular (IA) injections. Several actions should be taken prior to injecting Synolis VA: -Synolis VA gel should be at room temperature at the moment of the IA injection

- -The injection site must be carefully disinfected -Appropriate size of the needle must be used
- (recommendation for injection in the knee: 18 to 21 G, 2")
- -The needle must be firmly attached to the luer lock collar of the syringe
- -Inject accurately into the joint cavity only

When treating osteoarthritis of the knee, one injection of SYNOLIS VA could be sufficient for patient with low disease severity (Kellgren-Lawrence grades I/II &/or mild to moderate pain level). A regimen of 3 injections performed one week apart is particularly recommended for patients with high disease severity (Kellgren-Lawrence grades III/IV & severe to very severe pain level).

The safety and treatment efficacy of SYNOLIS VA beyond 6 months have not been established.

Contra-indications

Synolis VA must not be:

- -used in patients with known hypersensitivity to sodium hyaluronate and/or sorbitol
- -injected into an infected joint
- -injected in patients with a skin disorder or an infection at or near the site of the injection
- -injected intravascularly
- -injected in pregnant or breast feeding women
- -injected in young people under the age of 18 years -injected into a joint of a limb presenting venal or lymphatic
- stasis
- -injected into an infected or seriously inflammed joint

Precautions for use

Before treatment the patient must be informed about the device, its contra-indications and possible side effects.

- In the absence of available clinical data on tolerance and efficacy of Synolis VA in patients with antecedents or active auto-immune disease, or patients with an abnormal physiological condition, the physician must decide whether to inject Synolis VA on a case-by-case basis depending on the nature of the disease as well as the associated concomitant treatments. It is also recommended to carefully monitor these patients after injection.
- -Check the integrity of the inner packaging prior to use and check the expiry date. Do not use the product if the expiry date has lapsed or if the packaging has been opened or damaged.
- -Do not transfer Synolis VA inot another container and do not add other ingredients to the product.
- -The IA injection should be performed carefully in order to avoid injecting outside the intra-articular cavity or into the synovial membrane. Viscoelastic gels injected in the persynovial area can be painful due to compression on the surrounding tissues.
- -In case of significant joint effusion, the physician must decide whether to inject Synolis VA on a case by case basis. Effusion must be aspirated before injecting Synolis VA.
- -Synolis VA is a single-use product, thus it should not be used for several patients and/or different sessions. The product must not be resterilised. Reuse of single-use products may cause infections as the sterility is void.
- -Only the gel is sterile, but not the outside of the syringe.
- -Synolis VA must be administered under strict aseptic procedures.
- -The patient is advised to avoid intense physical activity for at least 48 hours after the injection.
- -Product must be stored under recommended storage conditions

Drug interactions

There is a known incompatibility between sodium hyaluronate and quaternary ammonium salls such as benzalkonium chloride. Therefore, Synolis VA must never come into contact with such products (e.g.: certain disinfectants), nor with medical or surgical equipment treated with these types of products. To date, no data is available on the compatibility of Synolis VA with other products for intra-articular use.

Side effects

Intra-articular injections of Synolis VA may cause temporary pain, oedema and/or effusion. These reactions usually clear up in a few days. If these symptoms persist for over a week, or if any other side effects occur, the patient must inform his or her physician. The physician will be able to provide appropriate treatment for these effects.

Storage

Store between 2 and 25°C. Protect from light and extreme cold.

<u>l</u>	Manufactured by.
LOT	Batch Number .
STERILE	Product Sterile. Sterilized using steam.
<u> </u>	Refer to the instructions for use.
2	Single-use product. Do not reuse.
Ω	Use by date.
<u>^</u>	Caution.
+2°C 1 +25°C	Storage temperature : 2°C - 25°C.
*	Keep away from sunlight

CE marking in compliance with the directive 93/42/EEC relating to medical devices.

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APTISSEN S.A. Chemin du champ des filles 36