SYNOLIS VA

HYALURONIC ACID 80mg / SORBITOL 160mg

 APTISSEN S.A.

Chemin du champ des filles 36
1228 Plan-Les-Ouates
Geneva/Switzerland

Tel: +41(0) 848 949 848
mail@aptissen.com
www.aptissen.com

Information as of 06/2021
Updated documentation is available
from Aptissen SA. in Switzerland.

CE
1639

PSL019-2

GB Instructions for use

03

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06

Instructions for use of SYNOLIS VA

GB

Description

Synolis VA 80/160 is a viscoelastic, sterile, apyrogenic, isotonic, buffered, 2% solution of sodium hyaluronate. The sodium hyaluronate used in Synolis VA 80/160 is obtained from bacterial fermentation and presents a high mean molecular weight (MW) of 2 MDaltons. Synolis VA 80/160 has a neutral pH of 6.8 - 7.4 similar to the synovial fluid.

High concentration and MW of sodium hyaluronate combined with a polyol (sorbitol) that limits its degradation confer the ability of this viscoelastic solution to restore joint lubrication and shock-absorbing properties, similar to healthy synovial fluid. Synolis VA 80/160 functions by restoring physiological and viscoelastic properties of the synovial fluid which has been lost progressively during the osteoarthritis (OA) development. Therefore Synolis VA 80/160 reduce local pain and discomfort caused by symptomatic OA and improves mobility of the synovial joints.

Introduction

Synolis VA 80/160 is available as a pre-filled 5ml glass syringe containing 4 ml of visco-analgesic gel. Synolis VA 80/160 is available in 1-syringe box packaging. Each box includes traceability labels and 1 instruction leaflet.

Composition

For 1 ml:

Sodium hyaluronate	20 mg
Sorbitol	40 mg
Phosphate-buffer <i>qs</i>	1 ml

Synolis VA 80/160 syringe contains 4 ml of viscoelastic solution. Sterilisation method: moist heat.

Indications for use

Synolis VA 80/160 is indicated for treatment of symptomatic osteoarthritis (OA), in order to reduce pain and improve mobility following degenerative changes in the knee synovial joint. This treatment is intended for patients who have failed conservative nonpharmacologic therapy and simple analgesics and/or NSAIDs or who have intolerance to simple analgesics and/or NSAIDs.

Dosage and method of administration

When treating osteoarthritis of the knee, Synolis VA 80/160 is administered as a single intra-articular injection.

An additional injection of Synolis VA 80/160 may be performed when OA symptoms resume or for maintaining local pain management and joint function. However, treatment benefits are expected to last for a minimum of 6 months for responding patients.

The time period before repeating the treatment regimen depends also on physician's experience and/or severity of the affection.

Synolis VA 80/160 should be injected within the synovial cavity by a skilled medical doctor trained on Intra-Articular (IA) injections. Several actions should be taken prior to injecting Synolis VA 80/160:

- Any joint effusion must be aspirated before injecting the viscoelastic gel.
- Synolis VA 80/160 should be at room temperature at the time of IA injection.
- The injection site must be carefully disinfected.
- Appropriate size of needle must be used (recommendation: 18G to 21G).
- The needle must be firmly attached to the luer lock collar of the syringe.
- Inject accurately into the joint cavity only.

Contra-indications

Synolis VA 80/160 must not be:

- Injected in patients with known hypersensitivity to sodium hyaluronate and/or sorbitol preparation;
- Injected in pregnant or breast-feeding women;
- Injected in patients under the age of 18 years;
- Injected intravascularly;
- Injected in patients with a skin disorder or an infection at the site of the injection.

Precautions for use

Before treatment, patients must be informed about the device, its contra-indications and possible side effects. Do not use Synolis VA 80/160 for any indication other than symptomatic OA. In the absence of available clinical data on tolerance of the injection of Synolis VA 80/160 in patients with antecedents or an active auto-immune disease, the physician must decide whether to inject Synolis VA 80/160 on a case-by-case basis, depending on the nature of the concomitant diseases as well as the associated concomitant treatments. It is recommended to propose a prior test to these patients and not to inject if the disease is evolving. It is also recommended to carefully monitor these patients after injection. Check the integrity of the inner packaging prior to use and check the expiry date. Do not use the product if the expiry date has lapsed or if the packaging has been opened or damaged. Do not transfer Synolis VA 80/160 into another container and do not add other ingredients to the product.

The intra-articular injection should be performed carefully in order to avoid:

- Injecting into a joint of a limb presenting venal or lymphatic stasis;
- Injecting into an infected and inflamed joint;
- Injecting outside the intra-articular cavity or into the synovial membrane. Viscoelastic gels injected in the peri-synovial area can be painful by compressing other surrounding tissues;
- Injecting if there is significant joint effusion.

Synolis VA 80/160 solution must be administered under strict aseptic conditions. Synolis VA 80/160 is a single-use product, which must not be used for different patients and/or different sessions. It must not be re-sterilised.

Reuse of single-use products may cause infections as the sterility is void. Only the gel is sterile but not the outside of the syringe. Patients are advised to avoid any intense physical activity for at least 48 hours after the intra-articular injection. Discard the syringe in accordance with accepted medical practice and applicable national, local and institutional requirements.

Drug interactions

There is a known incompatibility between sodium hyaluronate and quaternary ammonium salts such as benzalkonium chloride. Synolis VA 80/160 should not come into contact with such products (this concerns certain disinfectants), nor with medical or surgical equipment treated with this type of products. To date, no data is available on the compatibility of Synolis VA 80/160 with other products for IA injection use.

Side effects

Possible side effects exist and must be described to the patient before treatment. Slight bleeding may occur during the injection, although it stops spontaneously as soon as the injection is completed. In occasional cases one or more of the following reactions may occur either immediately or as a delayed reaction. It can be temporary local pain, oedema, and/or joint effusion. These reactions usually heal within few days. If these symptoms persist for over a week, or if any other side effects occur, the patient must inform the doctor. The doctor may prescribe appropriate treatment for these undesirable effects. Other possible typical side effects of viscosupplement injections include, inflammation, redness, swelling, skin irritation, allergic and tissue reaction.

Storage

Store between 2°C and 25°C. Protect from light and extreme cold. Do not freeze the product.

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